


## Disinfection Guidelines



The following are approved for use with RIVANNA Accuro® systems:

### Ultrasound Gels

Aquagel  
 Aquasonic 100  
 Aquasonic Clear  
 Clear Image  
 Clear Wave  
 Scan  
 Ultraphonic

### Cleaners/Disinfectants

Clinell 2% Chlorhexidine Wipes  
 Transeptic Spray  
 Sono Wipes  
 P.I. Spray II  
 Protex Wipes  
 Sani-Cloth HB Wipes  
 Sani-Cloth Plus Wipes  
 Sani-Cloth Bleach

### High-Level Disinfectants

Sterris VPro Max2  
 Sterris VPro S2  
 Sterrad NX  
 Sterrad 100NX  
 Sterizone VP4

